


# Advancing Eyelid Surgery Training with 3D Printed Anatomical Models





With these lifelike 3D models, we can simulate rare and complex pathologies, offer ethical and reproducible training, and push surgical education to the next level.”



**Alexander Hechenberger**  
CEO  
Addion GmbH



## The Challenge

### Providing Standardized, Hands-On Training in Eyelid Surgery

At the University of Innsbruck, leaders in plastic and reconstructive surgery are addressing a critical gap in surgical education: the lack of standardized, realistic models for eyelid surgery training.

Traditional training has long relied on body donors' tissue or animal models, which present challenges ranging from limited availability to ethical and cultural constraints. These methods also lack the ability to consistently simulate complex or rare pathologies, leaving trainees underprepared for high-stakes procedures like blepharoplasty, ptosis correction, and tumor resections.

“Body donors are invaluable, but they can't provide standardized pathologies or realistic tissue behavior due to preservation,” explains Prof. Dr. Marko Konschake, director of the Institute of Clinical and Functional Anatomy. “We needed an additional tool that could bridge the gap between anatomical accuracy and surgical realism.”

Enter Addion GmbH, a pioneer in 3D printed medical devices, and their collaboration with Stratasys. Together with eyecre.at, the company specializing in the development and distribution of 3D tissues and organs, the team set out to revolutionize the hands-on training experience at the 2025 Lid Summit Tirol, one of Europe's premier workshops for aesthetic and reconstructive eyelid surgery.

## The Solution

### High-Fidelity 3D Printed Eyelid Models

To meet the needs of the Lid Summit Tirol, Addion developed a specialized eyelid model using Stratasys' advanced J850™ Digital Anatomy™ Printer. Designed to replicate delicate eyelid tissues, skin elasticity, and adipose layers, the model offers unprecedented realism for surgical rehearsal.

The 3-day course brought together 20 surgeons and trainees from across Europe for immersive practice on both body donors and 3D printed eyelid models. With 10 workstations running 8 hours daily, participants performed multiple procedures per day—enhancing their skills through repetition and exposure to simulated pathologies, including tumors and trauma defects.

“The ability to repeatedly practice on a standardized model with realistic tissue behavior is a game changer,” says Prof. Dr. Anton Schwabegger, director of the University Clinic for Plastic, Reconstructive and Aesthetic Surgery. “It improves safety, skill, and confidence before ever stepping into the operating room.”





## The Results

### Safer, Smarter Surgical Training

Feedback from the summit was resounding: the 3D printed models elevated both learning outcomes and learner confidence. Participants praised the tactile realism, lifelike coloration, and ability to simulate complex scenarios rarely encountered in body donor training.

“The models not only complement body donor practice—they add entirely new dimensions,” says Prof. Dr. Konschake. “From advanced flap techniques to retrobulbar anesthesia, trainees could rehearse procedures repeatedly, accelerating their mastery of critical skills.”

Surgeons emphasized that 3D models also address logistical challenges. Unlike body donors, which are limited in supply and require special storage and transport, 3D printed models are portable, scalable, and can be used repeatedly—making them ideal for international workshops and remote training environments. Additionally, body donors often present unpredictable variables, which can affect the consistency and quality of surgical training. In contrast, 3D models provide a standardized experience, allowing for more controlled and effective practice sessions.

“You can bring them anywhere, store them for months, and always have consistent quality ready for training,” says Addion GmbH’s founder and CEO Alexander Hechenberger.

## The Technology

### J850™ Digital Anatomy™ Printer

Stratasys’ J850 Digital Anatomy Printer was at the heart of this breakthrough. Its unique capability to blend multiple photopolymers produces models with precise mechanical, visual, and tactile properties. From soft tissue to bone-like structures, it delivers lifelike simulation unmatched by conventional materials.



For the eyelid models, Addion engineers worked closely with Stratasys to fine-tune material properties, capturing subtle details like skin pliability, muscle tension, and fat pad behavior. The result: a model that behaves strikingly like living tissue when manipulated with surgical tools.

#### Key Anatomical Features to Simulate

- Skin pliability
- Muscle tension
- Fat pad behavior
- Pathological variations including tumors



3D printing is reshaping how we prepare surgeons for the real world. With each innovation, we take another step toward safer, more precise patient care.”

**Prof. Marko Konschake**  
Institute of Clinical and Functional  
Anatomy, Medical University of  
Innsbruck (MUI), Innsbruck, Austria



# The Impact

## Redefining the Future of Medical Training

The success of the Lid Summit Tirol points toward a new future for surgical education—one where 3D printed models serve as essential complements to body donor and live-patient training.

“We’re not stopping until every clinic has access to a 3D printing department,” says Hechenberger. “We envision a world where doctors routinely print patient-specific models for both training and procedural planning.”

For the University of Innsbruck, this collaboration marks just the beginning. Faculty plan to expand the use of 3D models into other surgical disciplines, including neurosurgery, otology, and trauma care.

“3D printing is reshaping how we prepare surgeons for the real world,” says Dr. Konschake. “With each innovation, we take another step toward safer, more precise patient care.”

### Materials and Methods

- Printer Used: Stratasys J850™ Digital Anatomy™ Printer
- Materials Used: Agilus30™ Colors and Digital Anatomy materials
- Collaborators: Addion.at GmbH, eyecre.at GmbH, MED-LAB GmbH, Medical University of Innsbruck (MUI), MUI - Institute of Clinical and Functional Anatomy, University Clinic for Ophthalmology and Optometry, and University Clinic for Plastic, Reconstructive and Aesthetic Surgery, and Cluster Life Sciences Tirol
- Application: Eyelid surgery training (blepharoplasty, ptosis correction, tumor resection)
- Training Context: 2025 Lid Summit Tirol (28 surgeons, 3-day course, 15 stations, 4 operations per day)

Learn how Stratasys Digital Anatomy™ solutions are transforming surgical training worldwide at [stratasys.com/en/industries-and-applications/3d-printing-industries/medical/](https://stratasys.com/en/industries-and-applications/3d-printing-industries/medical/).



**stratasys.com**

ISO 9001:2015  
Certified

Stratasys Headquarters  
5995 Opus Parkway,  
Minnetonka, MN 55343  
+1 800 801 6491 (US Toll Free)  
+1 952 937-3000 (Intl)  
+1 952 937-0070 (Fax)

1 Holtzman St., Science Park,  
PO Box 2496  
Rehovot 76124, Israel  
+972 74 745 4000  
+972 74 745 5000 (Fax)

### CASE STUDY

#### MEDICAL

©2025 Stratasys. All rights reserved. Stratasys, the Stratasys Signet logo, PolyJet, J850, and Vero are registered trademarks of Stratasys Inc. All other trademarks are the property of their respective owners, and Stratasys assumes no responsibility with regard to the selection, performance, or use of these non-Stratasys products. Product specifications subject to change without notice. CS\_PJ\_ME\_Eyelid Surgery\_0725a